

# Radixact™

A UNIQUELY SMART  
RADIATION THERAPY SYSTEM




Precise, innovative tumor treatments™

**Radixact™**

# FULLY INTEGRATED TREATMENT PLANNING AND DELIVERY

As the radiation oncology field advances, interoperability is becoming a critical need. The Radixact™ Treatment Delivery System provides a fully-integrated platform for uniquely smart treatment planning, data management and treatment delivery.



## Smart Treatment Applications

The Radixact™ Treatment Delivery System delivers megavoltage x-ray radiation using rotational, non-rotational, intensity-modulated radiation therapy, or non-intensity-modulated radiation therapy (three-dimensional conformal) treatment techniques, and using image-guided radiation therapy or non-image-guided radiation therapy workflows, in accordance with the physician-approved plan.

For the full range of radiation oncology patients

# A UNIQUELY SMART RADIATION THERAPY SYSTEM

## PLAN

Accuray Precision™ System

### KEY FEATURES

- Planning for multiple Radixact™ Systems
- Dose calculation for multiple systems on a single Accuray Precision™ workstation
- Treatment plan comparison
- AutoSegmentation™ options (head and neck, brain, prostate)
- Comprehensive contouring tools
- Remote review and approval of plans
- User-configurable dose volume display options
- Multi-modality image support
- PreciseART™ Adaptive Radiation Therapy Option
- PreciseRTX™ Retreatment Option

## SEE

iDMS™ Data Management System

### KEY FEATURES


- Remote review of daily imaging
- Digital signature capabilities
- Individual user privileges
- Patient record management for multiple Accuray Systems
- DICOM data import/export and Oncology Information System interfaces
- Data storage and management tools with storage vault option
- Disaster recovery features, including automatic backup and data mirroring

## TREAT

Radixact™ Treatment Delivery System


### KEY FEATURES

- 1,000 MU/min linac
- 10RPM gantry rotation for image acquisition
- 6MV system with integrated beamstop
- TomoHelical™ and TomoDirect™ modes
- TomoEDGE™ dose sculpting mode



# Radixact™

## A UNIQUELY SMART RADIATION THERAPY SYSTEM



Precise, innovative tumor treatments™

### UNITED STATES

#### Accuray Corporate Headquarters

1310 Chesapeake Terrace  
Sunnyvale, CA 94089  
USA

Tel: +1.408.716.4600  
Toll Free: 1.888.522.3740  
Fax: +1.408.716.4601  
Email: sales@accuray.com

#### Accuray Incorporated

1240 Deming Way  
Madison, WI 53717  
USA

Tel: +1 608.824.2800  
Fax: +1 608.824.2996

### ASIA

#### Accuray Japan K.K.

Shin Otemachi Building 7F  
2-2-1 Otemachi, Chiyoda-ku  
Tokyo 100-0004  
Japan

Tel: +81.3.6265.1526  
Fax: +81.3.3272.6166

#### Accuray Asia Ltd.

Suites 1702 - 1704, Tower 6  
The Gateway, Harbour City  
9 Canton Road, T.S.T.  
Hong Kong

Tel: +852.2247.8688  
Fax: +852.2175.5799

#### Accuray Accelerator Technology

#39 Huatai Road, Longtan  
Industrial Zone,  
Section 2 East, 3rd Ring Road,  
Chengdu, Sichuan 610051  
People's Republic of China

### EUROPE

#### Accuray International Sarl

Route de la Longeraie 9 (3rd floor)  
CH - 1110 Morges  
Switzerland

Tel: +41.21.545.9500  
Fax: +41.21.545.9501